




Easy to use
Quick to record
Convenient to link!



Horus ENT Scopes



Horus ENT Scope Set





Horus ENT Scope



PHYSICIAN



DIAGNOSTIC
STATION



SPECIALIST
NETWORK



TELEMEDICINE
& EMR



MOBILE
MEDICINE



REMOTE
MEDICINE

- Lightweight and easy to carry.
- Facilitating EMR and Telemedicine with Digital Image and Record.
- High resolution image and video to enhance doctor-patients communication.
- Good interoperability through WIFI, Micro-SD card, USB to transfer image for second opinion and remote consultation.
- Green features
Modular design (One control unit for different lens)/ Shortening the supply chain.


Horus ENT Scope Set

Otoscope	wide focus range 5mm-50mm
General lens	f-4mm, F/1.5
Resolution	Full HD true 1920x1080
Light Source	Natural White Light Emitting Diode
LCD Display	3.5" Full HD Color LED
Image Format	JPEG(Photograph) and H.264(Video)
File Transfer	Mini-USB 、AV output
File Storage	Micro SD Card (supports 2G to 32GB)
Power Source	Rechargeable Lithium Battery 3.7V
External Power	100-240 VAC, 50/60 Hz
Operating Time	3 hours
Weight	350g (approx.)



All in one
ENT scope

EOC 100



Product Name	Horus ENT Scope		
Green Features	One control unit for different lens/ Shortening the supply chain		
Origin	Taiwan		
Website	www.miis.com.tw	Service mail	sales@miis.com.tw
Manufacturer	Taiwan		
Customer Service	+886-3-579-8860		
Address	30076 1F, No. 7, R&D RdII, Hsinchu Science Park, Hsinchu, Taiwan		



GREEN TRADE PROJECT OFFICE
MINISTRY OF ECONOMIC AFFAIRS



GREEN MARKETING
& CONSULTING