



LiNA Port

Single use 5mm trocar system



Ordering information:

LiNA Port

REF.	Product: Dilating Tip	Packing
1005-06D1	5mm Cannula and Dilating Trocar - short	Single pack, sterile / Box of 5 sets
1005-06D2	5mm Cannula and Dilating Trocar - short	Dual pack*, sterile / Box of 5 sets
1005-09D1	5mm Cannula and Dilating Trocar - long	Single pack, sterile / Box of 5 sets
1005-09D2	5mm Cannula and Dilating Trocar - long	Dual pack*, sterile / Box of 5 sets
1005-06D1-GV	5mm Cannula with gas port and on/off valve and Dilating Trocar - short	Single pack, sterile / Box of 5 sets
1005-06D2-GV	5mm Cannula with gas port and on/off valve + one Cannula without gas port + Dilating Trocar - short	Dual pack*, sterile / Box of 5 sets
1005-09D1-GV	5mm Cannula with gas port and on/off valve and Dilating Trocar - long	Single pack, sterile / Box of 5 sets
1005-09D2-GV	5mm Cannula with gas port and on/off valve + one Cannula without gas port + Dilating Trocar - long	Dual pack*, sterile / Box of 5 sets
1005-06D1-G	5mm Cannula with gas port without on/off valve and Dilating Trocar - short	Single pack, sterile / Box of 5 sets
1005-06D2-G	5mm Cannula with gas port without on/off valve + one Cannula without gas port + Dilating Trocar - short	Dual pack*, sterile / Box of 5 sets
1005-09D1-G	5mm Cannula with gas port without on/off valve and Dilating Trocar - long	Single pack, sterile / Box of 5 sets
1005-09D2-G	5mm Cannula with gas port without on/off valve + one Cannula without gas port + Dilating Trocar - long	Dual pack*, sterile / Box of 5 sets

REF.	Product: Conical Tip	Packing
1007-1	5mm Cannula and Conical Trocar - short	Single pack, sterile / Box of 5 sets
1007-2	5mm Cannula and Conical Trocar - short	Dual pack*, sterile / Box of 5 sets
1007-2-2GP	5mm Cannula and Conical Trocar with 2 gas port - short	Dual pack*, sterile / Box of 5 sets
1009-C-1	5mm Cannula and Conical Trocar with gas port - short	Single pack, sterile / Box of 5 sets
1009-C-2GP	5mm Cannula and Conical Trocar with 2 gas port - short	Dual pack*, sterile / Box of 5 sets
2007-1	5mm Cannula and Conical Trocar - long	Single pack, sterile / Box of 5 sets
2007-2	5mm Cannula and Conical Trocar - long	Dual pack*, sterile / Box of 5 sets

REF.	Product: Pyramidal Tip	Packing
1005-1	5mm Cannula and Pyramidal Trocar - short	Single pack, sterile / Box of 5 sets
1005-2	5mm Cannula and Pyramidal Trocar - short	Dual pack*, sterile / Box of 5 sets
2005-1	5mm Cannula and Pyramidal Trocar - long	Single pack, sterile / Box of 5 sets
2005-2	5mm Cannula and Pyramidal Trocar - long	Dual pack*, sterile / Box of 5 sets
1008-1	5mm Cannula with gas port and on/off valve and Pyramidal Trocar - short	Single pack, sterile / Box of 5 sets
1008-2	5mm Cannula with gas port and on/off valve + one Cannula without gas port + Pyramidal Trocar - short	Dual pack*, sterile / Box of 5 sets
2008-1	5mm Cannula with gas port and on/off valve and Pyramidal Trocar - long	Single pack, sterile / Box of 5 sets
2008-2	5mm Cannula with gas port and on/off valve + one Cannula without gas port + Pyramidal Trocar - long	Dual pack*, sterile / Box of 5 sets
1009-1	5mm Cannula with gas port without on/off valve and Pyramidal Trocar - short	Single pack, sterile / Box of 5 sets
1009-2	5mm Cannula with gas port without on/off valve + one Cannula without gas port + Pyramidal Trocar - short	Dual pack*, sterile / Box of 5 sets
2009-1	5mm Cannula with gas port without on/off valve and Pyramidal Trocar - long	Single pack, sterile / Box of 5 sets
2009-2	5mm Cannula with gas port without on/off valve + one Cannula without gas port + Pyramidal Trocar - long	Dual pack*, sterile / Box of 5 sets

*2 x Cannula + 1 x Trocar + 1 Clear Sight Wipe



LiNA Medical ApS

• Formervangen 5
DK-2600 Glostrup
Denmark

• Tel.: +45 43 29 66 66
Fax: +45 43 29 66 99
info@lina-medical.com

• www.lina-medical.com
CVR: 10054974

LiNA Port

The LiNA Port is a single use trocar system fulfilling the demand for quality and efficiency at a competitive price.

safe

- **Silicone valve – airtight seal**
Ensuring minimum loss of pneumoperitoneum during surgery.
- **Dilating tip for delicate tissue**
The blunt dilating tip in the LiNA Port range supports non-traumatic entry.


functional

- **Threaded cannula**
To optimize efficient abdominal retention.

versatile

- **Choice of Dilating, Conical or Pyramidal tip**
Beside the dilating tip, the LiNA Port is available with conical or pyramidal tip for sharp and easy entry.
- **Two Lengths - Short and Long Ports**
Short 62mm and Long 84mm. Diameter inside cannula: Ø5,8mm allows easy insertion of all 5mm instruments.
- **Gas Ports**
The cannulas are available with either a simple gas port, a gas port with on/off valve or without gas port.
- **Single Pack and Dual Pack**
The LiNA Port is available in Single pack and Dual pack. The Dual Pack includes one trocar, two cannulas and one LiNA Clear Sight Wipe AntiFog Solution for camera optics.

Single Pack




Dual Pack with
one LiNA Clear
Sight Wipe



Airtight Silicone Double Seal Valve.

Designed for secure position of instrument, easy exchange of instruments and enables the use of instruments down to 3mm.



Non-traumatic Fascia Thread.
Designed for safe and firm anchoring of cannula.

Dilating, Conical or Pyramidal tip.
Supports the most used entry techniques.