


NEW PANTHER SMART STAPLER



Smart reload identification
Smart tissue recognition
Smart firing control



One touch firing
and opening
system

SMART HANDLE

Universal handle ready
to use with the most
extended portfolio of reloads

ELECTRONICALLY ASSISTED

Chip incorporated in the reloads
for the most precise firing



Innovative 45 degrees
Non-Stop articulation system

Universal handle with
better design and ergonomics
for easy 360 degrees knobbing



Technology P2G

- Double Gripping Technology
- Lateral diffusion of the fluids
- Better tissue compression

Curved tip Reloads
specially indicated for thin
tissue and vessels



The P2G technology
range guarantees
better stapling delivering
hemostasis and safety

NEW GENERATION ENDO STAPLING SYSTEM

Product Name	Product	Description		Color
Smart Stapler	PEAS	Powered Endoscopic Cutter Stapler with Non-Stop articulation system (Short Handle)		
	PEAM	Powered Endoscopic Cutter Stapler with Non-Stop articulation system (Standard Handle)		
	PEAL	Powered Endoscopic Cutter Stapler with Non-Stop articulation system (Extra long Handle)		
Endo Articulating Reloading	CADF30T	PANTHER Endo Articulating reloading unit 30 mm with P2G Technology and 2.5mm staple height (White)	Articulated reload	○ White
	CADF30N	PANTHER Endo Articulating reloading unit 30 mm with P2G Technology and 3.5mm staple height (Blue)	Articulated reload	● Blue
	CADF45T	PANTHER Endo Articulating reloading unit 45 mm with P2G Technology and 2.5mm staple height (White)	Articulated reload	○ White
	CADF45N	PANTHER Endo Articulating reloading unit 45 mm with P2G Technology and 3.5mm staple height (Blue)	Articulated reload	● Blue
	CADF-45S	PANTHER Endo Articulating reloading unit 45 mm with P2G Technology and 4.2mm staple height (Gold)	Articulated reload	● Gold
	CADF-45D	PANTHER Endo Articulating reloading unit 45 mm with P2G Technology and 4.8mm staple height (Green)	Articulated reload	● Green
	CADF-45R	PANTHER Endo Articulating reloading unit 45 mm with P2G Technology and 5.0mm staple height (Black)	Articulated reload	● Black
	CADF-60T	PANTHER Endo Articulating reloading unit 60 mm with P2G Technology and 2.5mm staple height (White)	Articulated reload	○ White
	CADF60N	PANTHER Endo Articulating reloading unit 60 mm with P2G Technology and 3.5mm staple height (Blue)	Articulated reload	● Blue
	CADF60S	PANTHER Endo Articulating reloading unit 60 mm with P2G Technology and 4.2mm staple height (Gold)	Articulated reload	● Gold
CADF60D	PANTHER Endo Articulating reloading unit 60 mm with P2G Technology and 4.8mm staple height (Green)	Articulated reload	● Green	
CADF-60R	PANTHER Endo Articulating reloading unit 60 mm with P2G Technology and 5.0mm staple height (Black)	Articulated reload	● Black	