




Electrodes for transurethral resection (TUR)



Application

Within the framework of a transurethral resection of the prostate and the bladder.
Also available as a **TURis** (Plasma) variant for use in physiological saline.


Compatible with Olympus and Storz resectoscopes.

Available in the following versions:

- Loop
- Ball 3 mm und 5 mm
- Knife
- Loop for vaporization

TUR electrodes are supplied sterile-packaged and can be used immediately.

TURis (Plasma) electrodes



Compatible with **Olympus** resectoscope

<i>Description</i>	<i>Art No.</i>	<i>PU.</i>
TURis (Plasma) versions		
Loop version	ME-ORS-S	10
Loop version, wide	ME-ORSW-S	10
Loop version, 30° optic	ME-ORS30-S	10
Loop version, 45° bended	ME-ORS45-S	10
Ball version, 3 mm	ME-ORK3-S	10
Knife version, 45° bended	ME-ORM45-S	10