

SOLUSCOPE TROLLEY & RACK


the next generation of storage
& transport solution



Soluscope

SOLUSCOPE TROLLEY &

Flexible



The Soluscope

Multiple use devices

ENHANCED MANAGEMENT OF YOUR ENDOSCOPES

The *Soluscope Trolley* and the *Soluscope Rack* are designed to optimise the endoscope turnover between examination and reprocessing rooms. They are compatible with *DSC8000* storage cabinet trays and *Soluscope* transport trays.

SOLUTION FOR ALL SITUATIONS

Beyond its transport function, the *Soluscope Trolley* is suitable to store up to 4 ready to use endoscopes for multiple exams or when examination room is distant from the reprocessing room. It can also be used as a stand at the patient's bedside thanks to its ideally positioned tabletop.

FOR ALL KINDS OF ENDOSCOPES

The *Soluscope Rack* temporarily store clean endoscopes or specific endoscopes not frequently used. Allowing easy handling of trays, it can also receive dirty endoscopes after exams while waiting for reprocessing.

RACK

Ergonomic



The Trolley & Rack offers a concept specially designed for endoscope handling and easy transport.
This solution improves your team's daily workflow in an efficient manner giving you comfort and peace of mind.

To make your life easier

MODULARITY

Soluscope Trolley is compact and adapted to small spaces. It can transport up to 4 endoscopes.

It includes three accessories you can use according to your need: a tabletop for surgical sheets, gloves or other medical instruments, a document holder to transport patient files, and a specific compartment to transport cleaning agent used to disinfect the tray.

Its large wheels made of anti-noise material make it very handy and easy to drive.

The trolley is also remarkable with its original and modern design.

10 ENDOSCOPES

The *Soluscope Rack* can store up to 10 trays. Its small footprint allows fitting in small spaces.

It's the ideal companion to *Soluscope AERs* and cabinet portfolio to extend your storage capacity.

EASY MAINTENANCE

The *Soluscope Trolley* and the *Soluscope Rack* are easy to clean with a surface detergent or disinfectant solutions.

transport & storage of endoscopes

Secure



y transport.



Safety above all

HIGH LEVEL PROTECTION

Once positioned in the tray and covered with its individual lid, the endoscope is protected from physical shocks or kinking and from exterior contamination during its transportation to the examination room. Medical staff are protected from any microbiological contamination.

Moreover, the *Soluscope Trolley* and the *Soluscope Rack* have specific sliding tracks to prevent the trays from falling down.

STABLE

The guiding wheels of the *Soluscope Trolley* are equipped with breaks providing better stability.



TECHNICAL SPECIFICATIONS - TROLLEY

Capacity	4 trays
Dimensions (W x D x H)	50 x 60 x 87 cm (with handle) 50 x 60 x 92 cm (with a tray on the top)
Weight empty	20 kg

TECHNICAL SPECIFICATIONS - RACK

Capacity	10 trays
Dimensions (W x D x H)	61 x 61 x 193 cm
Weight empty	40 kg

TECHNICAL SPECIFICATIONS - TRANSPORT TRAY

Capacity	1 endoscope
Dimensions (W x D x H)	55 x 48 x 15 cm
Weight empty / with lid	1.5 / 2.1 kg





SOLUSCOPE COMPLETE SOLUTION

Soluscope proposes a complete solution for flexible endoscope reprocessing.

It offers products for every step of the reprocessing cycle.

Each part is certified, easy to use and effective in order to reduce your daily load.



www.soluscope.com

Soluscope SAS : 100, rue du Fauge - Z.I. Les Paluds - 13400 Aubagne - France

☎ +33 (0)4 91 83 21 22 ☎ +33 (0)4 91 83 21 10 ✉ info@soluscope.com

Soluscope products are compliant with *European and international norms*, and are *CE marked*. *ISO 9001 version 2008 and ISO 13485 version 2003 certified*, *Soluscope* is continuously committed to satisfying its clients rigorous demands.

distributed by

