

A.4018L

Self·led·light®

sapi med
DIAGNOSTIC
LINE

A.4019L

A.4021L

A.4022L

A.4023L

A.4024L



CE

Disposable **Self-led-light** anoscopes can be used by the specialist both in clinic and in the operating theatre.

Made of transparent plastic material and supplied in individually sealed bags with led light built-in. Led lighting can be performed by pulling out the transparent strip. Guaranteed same light intensity for two hours usage.


Compatible, through adapter, to systems of image acquisition (cameras) **ALREADY AVAILABLE IN HOSPITALS** (for example operating theatres laparoscopic column).

ALL ANOSCOPES WITH BUILT-IN LED LIGHT ARE AVAILABLE IN THE STERILE VERSION.

TECHNICAL FEATURES

- Clear, transparent plastic proctoscopes
- Flash free
- Guaranteed atraumatic
- Not rechargeable led lighting
- Separable anoscope and led lighting for safe disposable
- Conforms fully with essential requirements of European Council Directive 93/42/EU
- Latex free
- Non sterile/Sterile

REF A.4018L / A.4018LST (sterile)




Proctoscope paediatric/diagnostic flute beaked with led lighting built-in

Tube length: 43/56mm

Tube diameter: 18mm

Quantity per box: 25 pcs.

REF A.4019L / A.4019LST (sterile)




Proctoscope paediatric/diagnostic trunk-end with led lighting built-in

Tube length: 54mm

Tube diameter: 18mm

Quantity per box: 25 pcs.

REF A.4021L / A.4021LST (sterile) **Blue Scope Midi**



Surgical/Examination Proctoscope trunk-end with led lighting built-in

Tube length: 70mm

Tube diameter: 21mm

Blu Scope midi has similar design characteristics to code A.4023L/A.4023LST.


Allows excellent control of the anal canal walls.

Quantity per box: 25 pcs.

Also available in the version without adhesive

REF A.4021.1L.

REF A.4022L / A.4022LST (sterile) **Blue Scope**



Surgical/Examination Proctoscope trunk-end with led lighting built-in

Tube length: 79mm

Tube diameter: 23mm

Blu Scope has similar design characteristics to code A.4023L/A.4023LST. Al-


lows excellent control of the anal canal walls.

Quantity per box: 25 pcs.

Also available in the version without adhesive

REF A.4022.1L.

REF A.4023L / A.4023LST (sterile)



Surgical/Examination Proctoscope flute beaked with led lighting built-in

Tube length: 55/88mm

Tube diameter: 23mm

Lute beak to perform biopsies, elastic ligatures and haemorrhoid sclerosis, in-

cision of intra-anal haematomas, excision of hypertrophic anal papillae, sutures


with haemostatic stitiche following polypectomy, I.R. coagulation etc.

Quantity per box: 25 pcs.

Also available in the version without adhesive

REF A.4023.1L.

REF A.4024L / A.4024LST (sterile)



Surgical/Examination Proctoscope trunk-end with led lighting built-in

Tube length: 65mm

Tube diameter: 23mm

They have similar design characteristics to code A.4023L/A.4023LST. Allows

excellent control of the anal canal walls.

Quantity per box: 25 pcs.

Also available in the version without adhesive

REF A.4024.1L.