

Orchid Spec - designed for quality  
Used by 90% of Dutch hospitals

[go to Shop](#)



reddot design award  
winner 2013



BioBased Plastic  
available

Faster  
Friendlier  
Smoother  
Stronger  
Safer

Listed as **Advanced Speculum**  
on National Health Service  
frameworks UK (**NHS**)

## Preferred by doctors

- smooth **single-handed operation** and locking
- front handle design angled backwards providing **greater freedom of movement** & **greater range of access**
- **hands free** after placement
- inwardly folded edges **facilitate optimal Cervix orientation**
- **no preheating** required
- **unbreakable** plastic
- excellent product strength through **optimised design**
- up to **50% cost reduction**

## Preferred by patients

- **smooth rounded edges** ensuring patient comfort
- gap design **prevents “pinching”**
- single use - **guaranteed clean**
- patient **temperature friendly**
- **clinical** white look



## Why white

Research by leading women's health specialists has demonstrated specific advantages.

**Misdiagnosis reduced** - (Diagnosing through a 'clear' Speculum regularly leads to a distorted view due to local tissue depression and discolouration.)

White surface design vastly **improves light reflection and visibility**  
Patients perceive white as **high quality, hygienic** and **friendly**

Orchid Spec 'clear' provides advantages in situations such as inspecting sutures or bleeding.

## Standard edition - patient friendly sizes

Available in Clean White & Crystal Clear plastic


[go to Shop](#)




Virgo



Medium



Wide



Long



BioBased Plastic  
available

## Special SX edition

### SX - with smoke extraction

- Facilitates precise horizontal manipulation of instruments due to the slim line smoke channel
- 360 access around the Cervix due to inward folded edges at the tip of the blades positioning and lifting the Cervix
- Available in 3 standard sizes & 2 Open-Sided



Medium SX

### Ideally suited for

- Colposcopy
- LEEP surgery
- Laser surgery

[go to Shop](#)



Wide SX




Long SX

## Special Open-Sided edition

### Orchid Open - Open on one side - New

- Instruments or catheter easily inserted or removed through the open side
- generous spacious access
- Available in sizes Medium & Wide



Open

[go to Shop](#)




winner 2016  
reddot design award

### Ideally suited for

- Hysteroscopy
- HSG / THL
- Endometrial Ablation
- IUD insertion

### Open SX - smoke extraction

- Available in sizes Medium SX Wide SX



Open SX