








Foreign Body Management

Product overview

Foreign Body


OCTOPUS Foreign body retrieval net



	Product Code	Order No.	Characteristics	Volume Net	Length	ø Working Channel	Unit
	FOR1-N1-20-180	502477	Hexagonal, framed, Paediatrics	3 cm ³	180 cm	2,0 mm	5
	FOR1-N1-30-220	502479	Oval, framed	6 cm ³	220 cm	2,8 mm	5
	FOR1-N1-40-220	502480	Reservoir function, framed	8 cm ³	220 cm	2,8 mm	5
	FOR1-N2-30-220	502481	Kite, framed	6 cm ³	220 cm	2,8 mm	5
	FOR1-N3-30-220	502482	Kite, framed, rotatable	6 cm ³	220 cm	2,8 mm	5

Foreign body protection cap





	Product Code	Order No.	Characteristics	Unit
	WS0420000	500944	Foreign body protection cap	10

Foreign Body


Foreign body forceps



	Product Code	Order No.	Characteristics	Opening width	ø Tube	Length	ø Working Channel	Unit
	FOR1-A1-23-230	502185	Alligator, metal spring, coated	9 mm	2,3 mm	230 cm	2,8 mm	10
	FOR1-A2-23-230	502284	Hybrid jaw configuration alligator/rat tooth, metal spring, coated	9 mm	2,3 mm	230 cm	2,8 mm	10
	FOR1-A2-23-260	502547	Hybrid jaw configuration alligator/rat tooth, metal spring, coated Balloon Enteroscopy	8 mm	2,3 mm	260 cm	2,8 mm	10
	FOR1-P1-18-230	502548	Pelican, coated, Paediatrics	9 mm	1,8 mm	230 cm	2,0 mm	10
	FOR1-P1-23-230	502546	Pelican, coated	15,5 mm	2,3 mm	230 cm	2,8 mm	10
	FOR1-H3-23-230	501946	Three-arm grasping forceps, PTFE catheter	-	2,3 mm	230 cm	2,8 mm	5

KOALA Polyptrap



	Product Code	Order No.	Characteristics	Unit
	TRP1-S2	502499	2 chambers	25

contact@medwork.com
www.medwork.com

Tel.: +49 (0) 91 93 / 500 900
Fax: +49 (0) 91 93 / 500 999

medwork GmbH
Medworkring 1
91315 Höchstadt/Aisch
Germany